

CODE	TEST NAME	CPT4 CODE
LAB1209	Q FEVER AB Container: Red-top tube or gel-barrier tube Temperature: Maintain specimen at room temperature.	86638
LAB560	QUAD SCREEN MATERNAL, SERUM Container: Serum Gel Tube Temperature: Room temperature	86336 84702 82105 82677
LAB2399	QUANTIFERON GOLD Container: The QuantiFERON(R)-TB Gold Plus collection kit contains the instructions for the collection and handling of (one each): (1) gray-top (with white ring), uncoated (nil); (2) green cap with white ring, TB1 tube; (3) yellow cap with white ring, TB2 tube; (4) purple top with white ring, mitogen-coated. A high altitude kit is also available for locations between 3350 and 6150 feet; a cap with a yellow ring differentiates the tubes. Temperature: Maintain specimen at room temperature for four tube collection kit. Stability: Kit MUST be incubated within 16 hours of collection.	86480
LAB24658	QUANTIFERON TB GOLD (RESULTS ONLY) Container: The QuantiFERON(R)-TB Gold Plus collection kit contains the instructions for the collection and handling of (one each): (1) gray-top (with white ring), uncoated (nil); (2) green cap with white ring, TB1 tube; (3) yellow cap with white ring, TB2 tube; (4) purple top with white ring, mitogen-coated. A high altitude kit is also available for locations between 3350 and 6150 feet; a cap with a yellow ring differentiates the tubes. Temperature: 1. Maintain specimen at room temperature prior to incubation at 37 degrees C. Stability: Kit MUST be incubated within 16 hours of collection. 2. After incubation, maintain specimen at room temperature (17 degrees C to 27 degrees C). Do NOT centrifuge refrigerate, or ship tubes on ice.	86480
LAB3115	QUETIAPINE (SEROQUE) LEVEL, SERUM Container: Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended. Temperature: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.	80299
LAB23381	QUETIAPINE, UR Container: Urine transfer tube Temperature: Maintain specimen at room temperature up to three days. Specimens may also be refrigerated or frozen.	80342