

CODE	TEST NAME	CPT4 CODE
LAB24882	EBV AB PANEL 1, CSF Container: CSF (Cerebrospinal Fluid) Temperature: Refrigerate	86663 86665 x2 86664
LAB1234	ECHINOCOCCUS ANTIBODY, IGG Container: Red-top tube or gel-barrier tube Temperature: Refrigerate or FREEZE	86682
LAB16066	ECHOVIRUS AB Container: Sterile container Temperature: FREEZE	86658 x5
LAB20406	EGFR Container: Red-top tube, gel-barrier tube, OR green-top (lithium heparin) tube. Do NOT use oxalate, EDTA, or citrate plasma. Temperature: Maintain specimen at room temperature.[1]	82565
LAB794	EHRlichia ANTIBODY PANEL Container: Red-top tube or gel-barrier tube Temperature: Room temperature	86666 x4
LAB16	ELECTROLYTE PANEL Container: Gel-Barrier Tube Temperature: Room temperature	80051
LAB16_OREGON	ELECTROLYTE PANEL - LABCORP PORTLAND ONLY Container: 5.0 mL Gold Top Tube (SST); Also acceptable: 6.0 mL Red Top Tube (Plain); 3.0 mL Light Green Top Tube (Lithium Heparin) Pediatric Collection: 600 uL Light Green Top Tube Microtainer; Also acceptable: 600 uL Gold Top Tube Microtainer Temperature: Refrigerated (2-8°C)	80051
LAB565	ELECTROLYTES, URINE RANDOM Container: Plastic urine container, no preservative Temperature: Maintain specimen at room temperature.	84300 84133 82436
LAB3095	ELECTROLYTES, URINE, 24HR Container: Plastic urine container, no preservative Temperature: Maintain specimen at room temperature.	84300 84133 82436
LAB24816	ENDOMYSIAL AB IGA + T-TRANSGLUTIMINASE IGA Container: Gel Barrier Tube	86364 x2 86231
LAB774	ENDOMYSIAL AB, IGA Container: Red-top tube or gel-barrier tube Temperature: Refrigerate	86231
LAB2090	ENTAMOEBA HISTOLYTICA AB Container: Red-top tube or gel-barrier tube Temperature: Refrigerate	86753
LAB1313	ENTAMOEBA HISTOLYTICA ANTIGEN Container: Unpreserved Stool Transport Vial Temperature: FREEZE - Critical frozen within 48 hours Stability: Room temperature unacceptable Refrigerate up to 48 hours Frozen up to 1 week	87337
LAB862	ENTEROVIRUS AB PANEL Container: Gel-barrier tube Temperature: FREEZE	86658 x5
LAB1333	ENTEROVIRUS, RNA, NAAT Container: Sterile container (CSF, stool), universal transport medium (NP/throat swab, rectal swab) Temperature: Swabs and stools are stable at room temperature, refrigerated, or frozen for seven days. CSF is stable at room temperature or refrigerated for seven days and frozen for 90 days.	87498