



Online Test Menu

LABORATORY CATALOG

LABCORP
WASHINGTON/OREGON
PROVIDENCE AMBULATORY

CODE	TEST NAME	CPT4 CODE
LAB303	FACTOR II, ACTIVITY, PROTHROMBIN SPECIFIC Container: Blue-top (sodium citrate) tube Temperature: Freeze	85210
LAB308	FACTOR IX, ACTIVITY Container: Blue-to (sodium citrate) tube Temperature: FREEZE	85250
LAB2078	FACTOR V ACTIVITY Container: Lavender-top (EDTA) tube, blue-top (sodium citrate) tube, yellow-top (ACD) tube OR buccal swab kit Temperature: Maintain specimen at room temperature. Stability: Room temperature up to 8 days Refrigerated up to 8 days Frozen up to 2 years	81241
LAB346	FACTOR V LEIDEN MUTATION Container: Lavender-top (EDTA) tube, yellow-top (ACD) tube, OR LabCorp buccal swab kit Temperature: Maintain specimen at room temperature.	81241
LAB304	FACTOR V, ACTIVITY Container: Blue-top (sodium citrate) tube Temperature: FREEZE	85220
LAB305	FACTOR VII ACTIVITY Container: Blue-top (sodium citrate) tube Temperature: FREEZE	85230
LAB2390	FACTOR VIII, ACTIVITY Container: Blue-top (sodium citrate) tube Temperature: Freeze Stability: Room temperature up to 4 hours Acceptable freeze/thaw up to 4 cycles	85240
LAB306	FACTOR VIII, ACTIVITY Container: Blue-top (sodium citrate) tube Temperature: FREEZE	85240
LAB1111	FACTOR VIII, MULTIMETRIC ANALYSIS Container: Blue-top (sodium citrate) tube Temperature: FREEZE	85247
LAB335	FACTOR VIII, RISTOCETIN COFACTOR Container: Blue-top (sodium citrate) tube Temperature: FREEZE	85245
LAB758	FACTOR X, ACTIVITY Container: Blue-top (sodium citrate) tube Temperature: FREEZE	85260
LAB309	FACTOR XI, ACTIVITY Container: Blue-top (sodium citrate) tube Temperature: FREEZE	85270
LAB310	FACTOR XII, ACTIVITY Container: Blue-top (sodium citrate) tube Temperature: FREEZE	85280
LAB1113	FACTOR XIII, QUAL Container: Blue-top (sodium citrate) tube Temperature: Freeze	85291
LAB390	FAT, FECAL, QUAL Container: Plastic screw-cap vial Temperature: Refrigerate at 2(degrees)C to 8(degrees)C.	82705