

LABORATORY CATALOG

CODE	TEST NAME	CPT4 CODE
LAB24788	JAK2 EXONS 12-15 REFLEX (RESULTS ONLY) Container: Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature: Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81279
LAB2127	JAK2 MUTATION PANEL Container: Lavender-top (EDTA) tube, green-top (sodium heparin) tube, Tan-top (K2-EDTA) tube, OR Pink-top (K2-EDTA) tube Temperature: Ship specimen at room temperature. Specimen should arrive in the laboratory within 48 hours of collection. If specimen is to be stored prior to shipment, store at 2(degrees)C to 8(degrees)C. Indicate date and time of collection on the test request form.	81270
LAB24732	JAK2 MUTATION QUAL WITH REFLEX TO CALR/MPL Container: Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature: Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81270
LAB24787	JAK2 V617F MUTATION ANALYSIS, QUAL, W REFLEX TO JAK2 EXON 12-15 Container: Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature: Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81270
LAB24998	JAK2 V617F MUTATION ANALYSIS, QUANTITATIVE Container: Lavender-top (EDTA) tube, green-top (sodium heparin) tube, Tan-top (K2-EDTA) tube, OR Pink-top (K2-EDTA) tube Temperature: Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81270
LAB24999	JAK2 V617F WITH REFLEX TO CALR + EXON 12-15 + MPL Container: Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature: Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81270
LAB24206	JC VIRUS AB W/INDEX, QUAL, REFLEX INHIB ASSAY Container: Red-top Tube OR Lavender Tube Temperature: Room temperature	NOCPT
LAB24850	JC VIRUS DNA, PCR, QUALITATIVE, CSF Container: Sterile container for CSF Temperature: Room temperature or refrigerated for up to 7 days. Frozen up to 90 days.	87798
LAB1888	JC VIRUS DNA, QUALITATIVE Container: Lavender-top (EDTA) tube	87798
LAB770	JO 1 AB, IGG Container: Red-top tube or gel-barrier tube Temperature: Room temperature	86235