

## LABORATORY CATALOG

CODE	TEST NAME	CPT4 CODE
LAB816	N-TELOPEPTIDE, CROSS-LINKED, URINE Container: Transfer tube, no preservative Temperature: Refrigerate	82570 82523
LAB13255	NATURAL KILLER CELL Container: Lavender-top (EDTA) tube AND yellow-top (ACD-A) or (ACD-B) tube Temperature: Maintain specimens at room temperature. Do NOT freeze or refrigerate specimen.	86355 86357
LAB13524	NEISSERIA GONORRHOEAE, CONFIRM Container: Gen-Probe(R) APTIMA(R) swab or APTIMA(R) urine specimen transport; ThinPrep(R) or SurePath(TM) liquid cytology vials Temperature: Maintain specimen at room temperature or refrigerate (2 degrees C to 30 degrees C)	87591
LAB262	NEISSERIA GONORRHOEAE, NAAT Container: Gen-Probe(R) APTIMA(R) swab or APTIMA(R) urine specimen transport; ThinPrep(R) liquid cytology vial Temperature: Room temperature	87591
LAB23734	NEOPTERIN Container: Red-top tube, gel-barrier tube or lavender-top (EDTA) tube Temperature: Refrigerate. Protect from exposure to ultraviolet light and sunlight.	83520
LAB2193	NEUROMYELITIS OPTICA AB Container: Red-top tube OR gel-barrier tube Temperature: Refrigerate	86051
LAB13263	NEURON SPECIFIC ENOLASE AG Container: Red-top tube or gel-barrier tube Temperature: Refrigerate	86316
LAB23621	NEURONAL CELL AB Container: Red-top tube or gel-barrier tube Temperature: FREEZE Stability: Refrigerate up to 7 days Frozen	83520
LAB1795	NEWBORN METABOLIC SCREEN REPEAT Container: PKU card Temperature: Room temperature	NOCPT
LAB413	NICKEL, URINE, 24HR Container: Plastic urine container, no preservative Temperature: Maintain specimen at room temperature.	83885 82570
LAB2435	NICOTINE AND METABOLITES Container: Lavender-top (EDTA) tube, green-top (heparin) tube OR red-top tube Temperature: Room Temperture	80323
LAB736	NICOTINE AND METABOLITES, URINE, QUANT Container: Urine container, no preservatives Temperature: Room Temperature	80323
LAB414	NICOTINE, URINE, QUAL Container: Urine container Temperature: Room Temperature	80307
LAB24894	NMO/AQP4 FACS, SERUM Container: Gel-barrier tube Temperature: Refrigerate	86053
LAB24278	NMR LIPOPROFILE TEST (LP-IR VALUES) Container: NMR LipoTube (black-and-yellow-top tube) is the preferred container, plain red-top tube, lavender-top (EDTA-no gel) tube, OR green-top (heparin-no gel) tube. Temperature: Refrigerate all acceptable tube types as soon as possible after centrifugation and within 24 hours of collection. Keep refrigerated prior to shipment, and SHIP ON FROZEN COOL PACKS. Do NOT store at room temperature. Do NOT freeze the sample. Sample is stable refrigerated for six days.	80061 83704