



# Online Test Menu

## LABORATORY CATALOG

LABCORP  
WASHINGTON/OREGON  
PROVIDENCE AMBULATORY

CODE	TEST NAME	CPT4 CODE
LAB13275	OBSTETRICS PANEL (REF) Container: gray top vacutainer for urine culture, urinalysis transport (red/yellow), gel-barrier tube, and lavender (EDTA)	87340 86850 85025 86762 81001 86900 86592 87389 86901 87086
LAB2502	OCCULT BLOOD, FECAL IMMUNOCHEMICAL TEST (FIT) Container: Polymedco OC-Auto Sampling Bottle Temperature: Specimen in sampling bottle can be stored for 15 days at ambient temperature. degrees C).	82274
LAB696	OCCULT BLOOD, GASTRIC Container: Sterile leakproof container with screw-cap. Refer to the Microbiology Specimen Collection chart in the Microbiology Appendix. Temperature: Refrigerate	82271
LAB3109	OLANZAPINE (ZYPREXA) LEVEL, SERUM Container: Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended. Temperature: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.	80299
LAB740	OLIGOCLONAL BAND PANEL Container: Red-top tube or gel-barrier tube AND plastic (CSF) tube Temperature: Room temperature	83916
LAB24928	ONCOLOGY FISH 1P, 19Q Container: Green-top (sodium heparin) tube; pediatric Vacutainer(R) is optimal Temperature: Maintain specimen at room temperature.	NOCPT
LAB7078	OP EXAM, FORMALIN ONLY Container: Para Pak-Pink Temperature: Maintain specimen at room temperature.	87177
LAB25125	OPIATES CONFIRMATION, URINE (NON ORD) Container: Use plastic urine drug bottle and tamper-evident seal for forensic specimen. Collection kits are available by request from the laboratory. Temperature: Maintain specimen at room temperature. Stability: If arrival at lab will extend beyond 7 days then refrigerate.	80361
LAB24424	OPIATES, CONFIRM, URINE Container: Plastic urine container Temperature: Refrigerate	80361
LAB2381	ORGANIC ACID, URINE Container: Plastic urine container, no preservative Temperature: FREEZE	83919
LAB418	ORGANIC ACIDS, URINE Container: Plastic urine container Temperature: FREEZE	83918
LAB107	OSMOLALITY, SERUM Container: Red-top tube, gel-barrier tube or green-top (heparin) tube Temperature: Refrigerate	83930
LAB201	OSMOLALITY, STOOL Container: Small plastic container; do NOT use large paint can. Temperature: FREEZE	84999
LAB420	OSMOLALITY, URINE Container: Plastic urine container, no preservative Temperature: Refrigerate	83935
LAB13444	OSTEOCALCIN Container: Red-top tube or gel-barrier tube Temperature: FREEZE	83937